


# Novel neprilysin-derived BNP fragment in the circulation: evidence from a rat model

Feygina E.<sup>1,2</sup>, Artemieva M.<sup>2</sup>, Postnikov A.<sup>1,2</sup>, Tamm N.<sup>1,2</sup>, Bloschitsyna M.<sup>1,2</sup>, Medvedeva N.<sup>2</sup>, Katrukha A.<sup>1,2</sup>, and Semenov A.<sup>1,2</sup>  
1. HyTest Ltd. (Turku, Finland); 2. School of Biology, Moscow State University (Moscow, Russia)



## Background



**ENTRESTO™** IS A FIRST-IN-CLASS ARNI (ANGIOTENSIN RECEPTOR NEPRILYSIN INHIBITOR) DRUG FOR HEART FAILURE (HF) THERAPY, WHICH COMPRISES NEPRILYSIN (NEP) INHIBITOR SACUBITRIL.

ONE OF THE IMPORTANT NEP SUBSTRATES IS AN ACTIVE FORM OF B-TYPE NATRIURETIC PEPTIDE (BNP); THE AUGMENTATION OF THE ACTIVE BNP LEVEL DUE TO NEP INHIBITION IS CONSIDERED AS A POSSIBLE MECHANISM OF CARDIAC FUNCTION IMPROVEMENT BY ARNI.

WE HYPOTHEZIZED THAT PRODUCTS OF BNP PROTEOLYSIS BY NEP ARE PRESENT IN THE CIRCULATION AND THAT THEIR LEVEL MIGHT REFLECT THE NEP IMPACT ON THE ACTIVE BNP METABOLISM.

WE SUGGEST THAT BNP CLEAVAGE BY NEP AT 17-18 AAR RESULTS IN THE BNP RING STRUCTURE OPENING AND THE FORMATION OF THE **BNP-NEO17 FORM**.

### AIM OF THE STUDY


EXPLORE THE NEP-DEPENDENT **BNP-NEO17** GENERATION  
*IN VITRO* AND *IN VIVO*

## Methods

### IMMUNOASSAYS

#### BNP-NEO17 SANDWICH

**POLY-NEO17**  
POLYCLONAL RABBIT ANTIBODY, SPECIFIC TO NEO17-EPIPEPE



**50E1**  
MOUSE MONOCLONAL ANTIBODY (HYTEST), EPIPEPE 26-32

#### TOTAL BNP ASSAY


**SES (SINGLE EPIPEPE SANDWICH)**

- MAB 24C5 (EPIPEPE 11-17, HYTEST)
- MAB AB-BNP-2 (SPECIFIC TO THE IMMUNE COMPLEX OF BNP WITH MAB 24C5)




#### ANP COMPETITIVE ASSAY

**ANTI-ANP MAB 23/1**  
MOUSE MONOCLONAL ANTIBODY (BIORAD)



### IN VITRO BNP CLEAVAGE




- SYNTHETIC BNP 1-32 WAS INCUBATED WITH HUMAN AND RAT NEP EITHER WITH OR WITHOUT SACUBITRIL.
- SAMPLES WERE ANALYZED BY THE **BNP-NEO17 SANDWICH IMMUNOASSAY** AT 0, 0.5 AND 1 HOUR TIME POINTS.

**INCUBATION**

- 37°C
- 0.5h, 1h
- +/- SACUBITRIL


### IN VIVO BNP CLEAVAGE




- 8 WISTAR RATS WERE ANESTHETIZED AND THE RIGHT FEMORAL ARTERY AND VEIN CANNULATED.
- THE SAME RATS WERE USED FOR (-) AND (+) SACUBITRIL KINETICS.
- ANP WAS INJECTED TOGETHER WITH BNP AS A CONTROL FOR NEP INHIBITION AS A PREFERABLE NEP SUBSTRATE.
- SAMPLES WERE ANALYZED BY **BNP-NEO17, TOTAL BNP AND ANP IMMUNOASSAYS** AT 0, 1, 2, 4, 8 AND 20 MINUTES POST INJECTION

## Results

### BNP CLEAVAGE BY RAT AND HUMAN NEP *IN VITRO*



### CHANGE IN AVERAGE BNP AND ANP CONCENTRATIONS IN RAT PLASMA WITH AND WITHOUT NEP INHIBITION



## CONCLUSIONS

We have shown for the first time that **BNP-Neo17** is generated *in vitro* and *in vivo*, and that its formation is NEP-dependent. In the light of this, we suggest that **BNP-Neo17** concentration might reflect the inherent NEP activity of the organism and serve as a potential cardiac biomarker.