




Neuroscience

Antibodies and antigens



Introduction

Neurological disorders are ranked as the leading cause of disability-adjusted life years (DALYs) globally, with dementia being the third largest contributor of neurological DALYs (1). The World Health Organization (WHO) has endorsed a global action plan to react to the issue of dementia, which is an umbrella term for several diseases involving deteriorated cognitive functions and behavioral difficulties (2). Alzheimer's disease (AD) is the most common form of dementia, accounting for 60 - 70% of global cases (3). The WHO reported the number of people with dementia as being over 55 million worldwide in 2020. This figure is expected to reach 78 million in 2030 and 139 million in 2050 (3). In 2019, Alzheimer's disease and other forms of dementia ranked as the 7th leading cause of death (3).

Neurodegenerative diseases and traumatic brain injuries

Neurodegeneration is a highly complex process that is characterized by a progressive and irreversible loss of neurons from the specific brain and spinal cord regions (4). It is multifactorial and often associated with various molecular mechanisms, consequently leading to cognitive disability and dysfunctions. In neurodegenerative diseases, AD and Parkinson's disease (PD) are the most common ones. These predominantly affect the elderly population. The likelihood of developing a neurodegenerative disease rises dramatically with age.

Traumatic brain injury (TBI) and mild traumatic brain injury (mTBI) refers to injuries caused by an external force to the head. Previously, it was considered to cause a static neurological insult.




Figure 1. Neuroscience biomarkers and their locations.

Based on later studies, it has shown the potential to influence brain volume loss, which may set the stage for acute and also chronic neuroinflammatory effects, neuron damage, neurodegeneration, and dementia (5, 6). The synergistic effects of age and injury will accelerate the cognitive decline, executive dysfunction, resulting in the condition of dementia (6).

The relationship between neurodegenerative diseases and TBI shows critical health implications. The risk of chronic dementia and other neurodegenerative diseases is significantly elevated after TBI (6). Among the various neurological disorders, the biomarkers associated with the development of neurodegenerative diseases have often been relevant to TBI, such as for disease progression tracking.

Paving the way to enable early diagnosis

One of the critical challenges in terms of neurodegenerative diseases or neurological disorders is the early detection and correct identification of the causes of various neuro-logical conditions. The common diagnosis methods include cognitive tests, neuroimaging techniques such as magnetic resonance imaging (MRI) and positron emission tomography (PET), and detection of biomarkers from body fluids (such as cerebrospinal fluid, CSF, and blood) (6). The diagnosis of mTBIs and concussions is difficult as the imaging techniques fail to detect the case, while computed tomography (CT) may expose the patient to radiation and comes at a high cost.

Recent studies have shown the trend towards identifying novel biomarkers from samples of body fluids (4). However, diagnostic and prognostic values are still to be further discovered and fully established.

Hytest as a leading manufacturer of immunological reagents

The improvement in quality of life of older people is becoming a public health priority. As a leading reagent manufacturer, Hytest has been committed to developing high-quality immunological reagents for neurodegenerative diseases and other brain injuries since 1994. Hytest's neuroscience products include Neurofilament light (NfL), β -amyloid 42 (A β 42), Glial fibrillary acidic protein (GFAP), and S100 proteins, which are illustrated in Figure 1.


More information regarding the performance of the products and recommendations on capture-detection antibody pairs (when available) can be found on our website: www.hytest.fi

You may also contact our Sales Team directly by writing to: hytest@hytest.fi.

Abbreviations

AD	Alzheimer's disease
A β	β -amyloid, beta-amyloid
CLIA	Chemiluminescence immunoassay
CNS	Central nervous system
CSF	Cerebrospinal fluid
CT	Computed tomography
DALYs	Disability-adjusted life years
FDA	Food and Drug Administration
GFAP	Glial fibrillary acidic protein
Ifs	Intermediate filament
MAB	Monoclonal antibody
MRI	Magnetic resonance imaging
NfL	Neurofilament light
NSE	Neuron-specific enolase
PET	Positron emission tomography
PD	Parkinson's disease
TBI	Traumatic brain injury
mTBI	Mild traumatic brain injury
UCH-L1	Ubiquitin C-terminal hydrolase-L1
WHO	World Health Organization

Introduction	2	S100 proteins	12
Abbreviations.....	4	Additional products.....	14
Neurofilament light (NfL).....	6	Neuron-specific enolase (NSE).....	14
Beta-amyloid 1-42, human	8	Calmodulin	14
Glial fibrillary acidic protein (GFAP)	10	References.....	15



Neurofilament light (NfL)

CLINICAL UTILITY

- **Traumatic Brain Injury (TBI)**
- **Neurodegenerative diseases**

Neurofilaments are the main cytoskeletal components in neurons, which show their abundance in axons. They are bundles of 10 nm diameter filaments that are also classed as the intermediate filament (IF) proteins, which consist of a triplet subunit, including neurofilament light (NfL), neurofilament middle (NfM), and neurofilament heavy (NfH). Furthermore, there is a fourth subunit; either α -internexin in the central nervous system or peripherin in the peripheral nervous system, respectively (see Figure 2) (7).

NfL as a diagnostic marker of neuronal damage

Forming the core backbone of the neurofilament, NfL is essential to the structure of the axons and it has been found to have prognostic value. In recent years, it has emerged as a biomarker for neuronal damage and other neurologic conditions. Elevated NfL concentrations have been widely reported not only in CSF but also in blood under the condition of neurologic disorders, including multiple sclerosis, TBI, and stroke (8, 9).

Monoclonal antibodies specific to human NfL

Hytest provides several monoclonal antibodies (MAbs) that are specific to human NfL and that are derived from three different animal species and are designed to be used in NfL-specific immunoassays. Our antibodies do not have cross-reactivity to type III IF proteins (GFAP, vimentin, desmin, and peripherin)

and several other neuronal proteins. Pairs of MAbs effectively recognize both recombinant and endogenous NfL, and they may be used for a variety of immunoassays, such as direct EIA, indirect EIA, and sandwich-type immunoassays.

For more information about recommended pairs, please visit www.hytest.fi.





Figure 2.
Subunits of neurofilaments.

Detection of NfL in clinical samples

A prototype assay based on Hytest’s in-house pair recommendation was developed and used for the immunodetection of NfL in CSF samples from several patients with different diagnoses. NfL levels have been grouped into low, moderate, or high, as presented in Figure 3.

Figure 3.

NfL immunodetection in CSF samples. Prototype 3-step in-house sandwich-immunoassay NF79-NF71. The adjacent table lists the diagnosis of each patient. CSF samples were diluted 4-5-fold prior to measurement.



MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4NF3	Monoclonal anti-neurofilament light (NfL), human	Enzyme immunoassays

* Several MABs available under one catalogue number. Please see www.hytest.fi.

Beta-amyloid 1-42, human

CLINICAL UTILITY

- **Alzheimer's disease (AD)**

Alzheimer's disease (AD) is a complex progressive neurodegenerative disease and it is the most common cause of dementia. The three stages of AD include preclinical Alzheimer's, mild cognitive impairment (MCI), and dementia due to Alzheimer's (10). The importance of early diagnosis has been raised in both research and clinical societies before irreversible brain damage has occurred. The accumulation of the extracellular beta-amyloid protein (which is also referred to as A β plaques) is one fundamental neuropathological hallmark of AD. To date, A β is considered as being one of the most established biomarkers for the progression of AD.

Clinical practice of β -amyloid 1-42

Appropriate use criteria for CSF β -amyloid testing A β 42, in ratio with A β 40, was published in 2018 (11). Besides CSF, blood-based β -amyloid has also shown promise in the diagnosis of AD in recent years (12). Among the various A β isoforms found in CSF, A β 42 is considered as one of the key pathological

biomarkers in terms of assessment for AD diagnosis. A β 42 forms aggregation into fibrils and extracellular plaques in the brain, resulting in a lower level of A β 42 in the CSF. Therefore, A β 42 provides predictive value for AD, which can be used in both the prodromal and dementia stage of AD.

Monoclonal antibodies specific to β -amyloid 1-42

At Hytest, we provide well-characterized human β -amyloid-specific mouse MAbs for the detection of A β 42 in human CSF. These antibodies were developed against synthetic peptides that correspond to fragments of the A β 42 sequence. Our prototype assay immunoassays BAM7cc-BAM113cc and BAM7cc-BAM120cc correlate well with the commercially available INNOTEST β AMYLOID (142) assay, which is shown in Figure 4 on next page.




Figure 4. Correlation studies between the immunoassays BAM7cc-BAM113cc (A) and BAM7cc-BAM120cc (B) and INNOTEST β -AMYLOID(1-42) assay. The correlation coefficients (Pearson) between the assays and the INNOTEST assay are provided in the picture.

MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4BA3	Monoclonal mouse anti-beta-amyloid, human	Enzyme immunoassays

* Several MAbs available under one catalogue number. Please see www.hytest.fi.

Glial fibrillary acidic protein (GFAP)

CLINICAL UTILITY

- Traumatic Brain Injury (TBI)

Glial fibrillary acidic protein (GFAP) is an astrocyte structural protein. It is a member of the type III intermediate filament that provides support and strength to cells and sustains the cell shape. GFAP proteins appear as bundled fibers and form the main intermediate filament found in specific brain cells, which are called astrocytes. Astrocytes, which are also known as astroglia, are star-shaped cells that represent approximately 30 - 40% of the cells in the central nervous system (CNS) (13). They are tightly involved in the blood-brain barrier and they have direct interactions with various cells in the nervous system.

Clinical use of GFAP

In 2018, GFAP, together with Ubiquitin C-terminal hydrolase-L1 (UCH-L1), received approval as a biomarker for concussion (mTBI) in adults by the U.S. Food and Drug Administration (FDA) (14). This represented a great milestone in the blood biomarker development for neurological disorders. Furthermore, it can be detected in serum already within the first few hours after a head injury (15).

An increasing number of studies have indicated that GFAP might be a useful biomarker for the differentiation between a hemorrhagic and an ischemic stroke. These studies have shown that GFAP increases in the case of a hemorrhagic stroke within two hours after stroke onset, with peaking taking place between 6 and 12 hours after stroke onset. Instead, in the case of an ischemic stroke, the GFAP levels in blood increase at a later time point (16).

Monoclonal antibodies specific to GFAP

Hyttest offers several well-characterized MAbs that are specific to GFAP and which may be used for the quantification of GFAP in serum, plasma, or CSF. We have recommendations on the capture-detection pairs using a sandwich chemiluminescence immunoassay (CLIA). For more information, please visit www.hyttest.fi.

Moreover, we have tested GFAP in plasma samples that were obtained from patients with either a hemorrhagic (N=5) or an ischemic (N=5) stroke using our recommended MAb pairs (GFAP83cc–GFAP81cc). All of the samples were taken within the first 12 hours following the injury. The prototype assay only detected GFAP in the plasma samples from patients who had suffered a hemorrhagic stroke. This is in line with the results from other studies and it suggests that GFAP can be used for discriminating these two types of strokes. In ischemic strokes, the level of GFAP should only increase at a later stage.




Figure 5. GFAP was measured in plasma samples from hemorrhagic and ischemic stroke patients by using the GFAP83cc–GFAP81cc prototype assay.

Recombinant GFAP

Hyttest recombinant GFAP is suitable to be used as a standard or calibrator in immunoassays. The antigen consists of amino acids 60-383 of human GFAP and it is expressed in *E. coli*. Purity of the antigen is over 90%. SDS-PAGE of recombinant GFAP reveals that it migrates as one band with apparent molecular weight of 34 kDa (see Figure 6).

It should be noted that GFAP is a fibrillar protein and prone to polymerization. Also, the recombinant GFAP tends to form dimers and thus the purified protein preparation likely always contains some amount of dimeric GFAP. Long-term incubation at positive temperatures increases percentage of dimeric form, and can, eventually, lead to further aggregation.




Figure 6. SDS-PAGE of recombinant GFAP fragment under reducing conditions in a gradient gel (10-20%). 2 µg of purified protein was loaded on the gel.

Recombinant GFAP tolerates freeze-thaw cycles

We tested the immunochemical stability of recombinant GFAP after repeated freeze-thaw cycles. The immunoreactivity, as measured using GFAP83cc–GFAP81cc prototype assay, did not change significantly after the protein has been subjected to ten freeze-thaw cycles (see Figure 7).




Figure 7. Freeze-thaw stability of recombinant GFAP. Aliquot of recombinant GFAP was subjected to ten freeze-thaw cycles. Immunoreactivity after first, third, fifth and tenth cycles was measured using a sandwich immunoassay (GFAP83cc–GFAP81cc).

MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4G25*	Monoclonal mouse anti-human glial fibrillary acidic protein (GFAP)	Enzyme immunoassays Western blotting Immunohistochemistry

* Several MAbs available under one catalogue number. Please see www.hyttest.fi.

ANTIGENS

Cat.#	Product	Source	Purity
8G47	Glial fibrillary acidic protein (GFAP), human, recombinant	Recombinant	>90%

S100 proteins

CLINICAL UTILITY

- Severity of brain injury and neuronal damage
- Traumatic and ischemic brain injury
- Diagnosis and prognosis of malignant melanoma

S100 proteins constitute a family of approximately 20 calcium-binding proteins. These small proteins (10-12 kDa) have 20-50% homology of amino acid sequences but differ by origin and functions and may serve as markers of different pathological states. In brain tissue, S100 proteins are represented mainly by the S100BB homodimer and S100A1B heterodimer of approximately 21 kDa. They are synthesized in astroglial cells and can be used as sensitive and reliable markers for central nervous system damage (17).

S100 as a diagnostic marker of brain injury

Due to their abundance in the astrocytes of the brain, S100 proteins are important neuro biomarkers, such as for glial injury and blood-brain barrier damage (18, 4). They have shown elevated results over the first few days after severe TBI (16). In

the Scandinavian guidelines for head injury management (19), S100B has been recommended as an option to rule out patients after mild head injuries without the need for CT. Furthermore, their involvement in relation to AD patients has also been explored (20). It plays an important role in CNS development and recovery after injury. However, their direct role in the disease mechanism is still to be fully established.

Human brain S100 proteins

Hyttest's S100 proteins are purified from human brain tissue by several chromatographic methods, including gel filtration and ion exchange chromatography. After native gel electrophoresis by Ornstein-Davis, the protein is presented by two bands that correspond to A1B and BB forms (Figure 8).




Figure 8.
Native gel electrophoresis of S100 proteins (by Ornstein – Davis).

Antigen loaded:

Lane 1: 2 µg

Lane 2: 5 µg

Gel staining: Coomassie brilliant blue R-250



MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4S37	Monoclonal mouse anti-human S100 proteins	Enzyme immunoassays Western blotting

* Several MABs available under one catalogue number. Please see www.hytest.fi.

ANTIGENS

Cat.#	Product	Source	Purity
8S9h	S100BB homodimer and S100A1B heterodimer, human	Human brain	>95%
8S9b	S100BB homodimer and S100A1B heterodimer, bovine	Bovine brain	>95%
8S9-2h	S100BB homodimer, human	Human brain	>95%
8S9-2b	S100BB homodimer, bovine	Bovine brain	>95%

Additional products

Neuron-specific enolase (NSE)

Neuron-specific enolase (NSE) is a dimeric enzyme that is located primarily in the neuronal cytoplasm. NSE is the $\gamma\gamma$ homodimer isoform of the enzyme enolase. This isoform is found in neurons and cells of neuroendocrine origin as well as platelets and erythrocytes. The concentration of NSE has been shown to be elevated in serum after head injury in correlation

with the degree of cell damage in CNS (21). It is a biomarker for the prognostication of neurological outcomes after cardiac arrest and has also been suggested to be a biomarker in TBI (22).

MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4N6*	Monoclonal mouse anti-neuron specific enolase (NSE)	Enzyme immunoassay Western blotting Immunohistochemistry

* Several MAbs available under one catalogue number. Please see www.hytest.fi.

ANTIGENS

Cat.#	Product	Source	Purity
8NS3	Neuron-specific Enolase (NSE)	Human brain	>95%

Calmodulin

Calmodulin (CaM) is a small calcium-binding protein that plays an important role in the calcium signal transduction. It is a primary sensor, which, once bound to Ca^{2+} , undergoes a conformational change, allowing it to regulate CaM-binding proteins (23). Both CaM and its binding proteins

are involved in multiple events in amyloid pathways and the tangle formation, which is a known factor for AD. CaM plays an important role in normal neuronal function and in many aspects of neurodegeneration (24).

ANTIGENS

Cat.#	Product	Source	Purity
8C10b	Calmodulin, bovine	Bovine brain	>95%
8C10h	Calmodulin, human	Human brain	>95%

REFERENCES

1. World Health Organization. Follow-up to the political declaration of the third high-level meeting of the General Assembly on the prevention and control of non-communicable diseases. Seventy-fifth World Health Assembly, 2022, URL: https://apps.who.int/gb/ebwha/pdf_files/WHA75/A75_10Add4-en.pdf
2. World Health Organization. Global action plan on the public health response to dementia 2017 – 2025. WHO publications, 2017, URL: <https://www.who.int/publications/i/item/global-action-plan-on-the-public-health-response-to-dementia-2017---2025>
3. World Health Organization. Dementia. WHO fact sheets, 2022, URL: <https://www.who.int/news-room/fact-sheets/detail/dementia>
4. **Balasa A. F. et al.** Body fluid biomarkers for Alzheimer's disease – an up-to-date overview. *Biomedicines*, 2020, 8, 10, pp 421.
5. **Dadas A. et al.** Biomarkers in traumatic brain injury TBI: a review. *Neuropsychiatric Disease and Treatment*, 2018, 14, pp 2989 – 3000.
6. **Bigler E. D.** Traumatic brain injury, neuroimaging, and neurodegeneration. *Frontiers in Human Neuroscience*, 2013.
7. **Yuan A. et al.** Neurofilaments and neurofilament proteins in health and disease. *Cold Spring Harb Perspect Biol*, 2017, 9, 4, pp 018309.
8. **Gaetani L. et al.** Neurofilament light chain as a biomarker in neurological disorders. *Journal of Neurology, Neurosurgery & Psychiatry*. 2019, 90, pp 870 - 881.
9. **Khalil M. et al.** Neurofilaments as biomarkers in neurological disorders. *Nature reviews neurology*, 2018, 14, pp 577 – 589.
10. Alzheimer's Association. Diagnostic criteria and guidelines, 2022. URL: https://www.alz.org/research/for_researchers/diagnostic-criteria-guidelines
11. **Shaw L. M. et al.** Appropriate use criteria for lumbar puncture and cerebrospinal fluid testing in the diagnosis of Alzheimer's disease. *Alzheimer's & Dementia*, 2018, 14, pp 1505 – 1521.
12. **Hansson O. et al.** The Alzheimer's Association appropriate use recommendations for blood biomarkers in Alzheimer's disease. *Alzheimer's Dementia*, 2022, 7, pp 1 – 18.
13. **Abdelhak A. et al.** Blood GFAP as an emerging biomarker in brain and spinal cord disorders. *Nature Reviews Neurology*, 2022, 18, pp 158 – 172.
14. U.S. Food and Drug Administration. FDA authorizes marketing of first blood test to aid in the evaluation of concussion in adults, 2018. URL: <https://www.fda.gov/news-events/press-announcements/fda-authorizes-marketing-first-blood-test-aid-evaluation-concussion-adults>.
15. **Papa L. et al.** Time course and diagnostic accuracy of glial and neuronal blood biomarkers GFAP and UCH-L1 in a large cohort of trauma patients with and without mild traumatic brain injury. *Jama Neurol*, 2016, 73, 5, pp 551 – 560.
16. **Foerch C. et al.** Diagnostic accuracy of plasma glial fibrillary acidic protein for differentiating intracerebral hemorrhage and cerebral ischemia in patients with symptoms of acute stroke. *Clinical Chemistry*, 2012, 58, pp 237–245.
17. **Nylen K. et al.** Serum levels of S100B, S100A1B and S100BB are all related to outcome after severe traumatic brain injury. *Acta Neurochir* 2008, 150, pp 221 – 227.
18. **Kapural M. et al.** Serum S-100beta as a possible marker of blood-brain barrier disruption. *Brain Res* 2002, 14, 940, pp 102 – 104.
19. **Uden J. et al.** Scandinavian guidelines for initial management of minimal, mild and moderate head injuries in adults: an evidence and consensus-based update. *BMC medicine*. 2013, 11, 50.
20. **Mrak R. E. and Griffin W. S.** The role of activated astrocytes and of the neurotrophic cytokine S100B in the pathogenesis of Alzheimer's disease. *Neurobiology of Aging*. 2001, 22, 6, pp 915 – 922.
21. **Ergün R. et al.** Prognostic value of serum neuron-specific enolase levels after head injury. *Neurol Res*. 1998, 20, 5, pp 418-420.
22. **Haque A. et al.** Neuron specific enolase: a promising therapeutic target in acute spinal cord injury. *Metab Brain Dis*. 2016, 31, 3, pp 487-495.
23. **O'Day D. and Myre M.** Calmodulin-binding domains in Alzheimer's disease proteins: extending the calcium hypothesis. *Biochemical and Biophysical Research Communications*. 2004, 320, 4, pp 1051 – 1054.
24. **O'Day D.** Calmodulin binding proteins and Alzheimer's disease: biomarkers, regulatory enzymes and receptors that are regulated by calmodulin. *Int J Mol Sci*. 2020 Oct 5;21(19):7344.



**Scientific
excellence**
for IVD

HYTEST LTD

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku, FINLAND
Tel. +358-2-512-0900 • E-mail: hytest@hytest.fi
HYTEST.FI

HYTEST NORTH AMERICA INC.

800 MacArthur Blvd, STE 100,
Mahwah, New Jersey 07430
Tel. 833-368-3489 • E-mail: hytestna@hytest.fi
HYTESTNA.COM

OUR BROCHURES

Product Catalog 2026 • Cardiac Markers •
Neuroscience • Inflammation •
Infectious Diseases • Metabolic
Syndrome • Thyroid Diseases • Veterinary